

Dopplers and Probes

Non-Directional Dopplers

Please be sure to indicate both style and frequency when ordering.
ex: pencil, skinny, adult flat, infant flat, & 8.1 MHz, 8.3 Mhz, 9.1 MHz, 3.9 MHz, etc.

Model 811-B (Speaker, rechargeable battery)	\$785.
Model 811-BTS (Speaker, rechargeable, auto off)	\$850.
Model 811-BL (Speaker, rchg., cautery suppressor)	\$850.
Model 810-A (Large speaker, meter, "D" batteries)	\$1,035.
Model 915-BL (2 & 8 MHz for vascular and air embolus) ...	\$1,285.

Probes

Please be sure to indicate both style and frequency when ordering.
ex: pencil, skinny, adult flat, infant flat, & 8.1 MHz, 8.3 Mhz, 9.1 MHz, 3.9 MHz, etc.

For optimal performance the probe frequency you order must match the frequency of the instrument you are ordering for.

8 and 9 MHz probes in standard 3/8 inch diameter pencil style, or infant and adult flat probes	\$150.
<i>Same as above but with 7' or 10' cable*</i>	\$165.
8 and 9 MHz 1/4 inch skinny pencil probes	\$170.
<i>Same as above but with 7' or 10' cable*</i>	\$190.
4 MHz pencil probes	\$170.
<i>Same as above but with 7' or 10' cable*</i>	\$190.
2 MHz obstetrical and chest probes	\$160.
<i>Same as above but with 7' or 10' cable*</i>	\$180.
PPG Sensor (Push-in single pin plug)	\$150.
PPG Sensor (5-pin twist-on plug & shielded cable)	\$160.
Microvascular Probe*	\$205.
Avian Probe*	\$180.

Avian Probe only available in 8.2, 8.5, 8.9, 9.2, 9.3, 9.5 & 9.7 MHz
* **Please Note:** *These items are special order and not returnable*

For our Computerized Flo-Lab Model

8 MHz Contoured Pencil Probe -8.3MHz	\$185.
Red PPG Sensor	\$150.
Yellow PPG Sensor	\$150.
3.9 Pencil 7' GEN2 Probe with Parks' Ultra Shield Technology**...	\$295.
8.1 Pencil 7' GEN2 Probe with Parks' Ultra Shield Technology**...	\$275.

**Gen2 Doppler probes can only be used with Flo-Labs manufactured after 8/1/2022 or on Flo-Labs that have been modified/upgraded.

Pocket Doppler

Model 841-A (8 MHz, speaker, auto off)	\$700.
--	--------

Directional Doppler without Recorder

Model 909 (4 & 8 MHz, 2 speakers, rechargeable battery)..	\$1,925.
---	----------

Directional Dopplers with Recorder

Model 1050-C (8 MHz Doppler and PPG)	\$4,995.
Model 1059-C (8 & 4 MHz Doppler, PPG, PVR and cuffs)	\$7,495.

* While Supplies Last *

Vascular and Obstetrical

Model 917 (2 & 8 MHz, rechargeable, speaker)	\$1,035.
(Comes standard in the #900 aluminum carry case with removable lid and room for accessories)	

Obstetrical Only

Model 614-B (Rechargeable battery, small speaker)	\$650.
---	--------

800 & 901 Carry Case

The 614-B, 811-B, 811-BTS, 811-BL and 915-BL are available with either the #800 or #901 aluminum carry case with removable lid.

#800 carry case - \$160.00 additional if ordered with the unit, \$180.00 if ordered later, without the unit.

#901 carry case - \$180.00 additional if ordered with the unit, \$200.00 if ordered later, without the unit.

DOPPLERS: The determining factor in choosing a probe type and operating frequency is the depth of the vessel. For best results from shallow vessels we recommend a probe with a frequency in the range 8 and 9 MHz (ex: 8.1 MHz, 8.2 MHz, 8.3 MHz, 9.1 MHz, 9.2 MHz etc.) Most people find it adequate for all their vascular needs.

For dual frequency directional Dopplers, we recommend a probe in the range of 4 (ex: 3.9 MHz, 4.1 MHz) and in the range of 8 or 9 MHz. The 4 MHz probe is larger ($\frac{1}{2}$ inch diameter) but there is considerable improvement in the detection of deeper arteries and veins, especially on obese or swollen limbs.

FLAT PROBES: Use flat probes for routine blood-pressure measurements on the arms. We have infant and adult sizes. You can use the probes interchangeably. Both are made in the range of 8 and 9 MHz frequencies.

AVIAN PROBE: Veterinary use, call for details.

NOTE: The probe frequency must match the frequency of the instrument. You cannot use a 2 or 4 MHz probe on an instrument tuned to 8 or 9 MHz. You can mix the probe style, but not the frequency.

SURGERY: The Model 811 series provides adequate service in surgery because you normally only want to judge blood flow quality by ear. If you want blood flow indicated on a meter for later comparison, we recommend our Model 810-A. If you also need a recording, we recommend the Model 1050-C. Generally, in surgery, use a probe with a frequency in the range of 8 or 9 MHz. Microvascular probe also available, call for details.

CARRY CASES: Our workhorse Doppler in pediatrics, surgery, ICU and the emergency room is the Model 811-B and its variations. Models 811-B, 811-BL, 811-BTS, 915-BL and 614-B come standard in a rugged aluminum box. These instruments are also available in an attractive aluminum case with convenient carry handle, removable lid and room for accessories.

#800 case (9.25" length x 7.25" width x 4.00" deep) large enough to hold a bottle of gel, probe and battery charger -- available in blue only.

#901 case (11.30" length x 8.40" width x 5.00" deep) has the same features as above, but with additional space for a manometer and a few cuffs -- available in off white only.

If you own an older Doppler in a plastic case you can order a metal case for it. When you order a metal case, be sure to specify what model it is for.

Manufactured in the U.S.A. by:

PARKS

800-547-6427 www.parksmed.com