

PARKS

Lower Arterial

Name: Sample. Patient A
Date of Birth: 11/24/1958

Age: 56

Patient ID: 123456
Gender: Female

Ref By: Smith. James MD
Examined By: Brown. Amy RVT, RDMS
Read By: Jones. William MD

Date: 12/27/2021 2:09 PM
Accession #: XO-A062220-03
Room:: OP

History

Smoke: Never
Diabetic: Yes
Gangrene: None
Prev Angioplasty: No

Hypertension: Yes
Claudication: None
Vascular Ulcers: None
Prev Vasc Surg: No

Hyperlipidemia: Yes
Rest Pain: None
Skin Color Chg: None

Additional history notes or patient comments may be entered here. History questions as well as the answers are custom configured by the user.





Segmental BP

	Right		Left
Brachial:	138	Index	142
Thigh:	162	1.14	169
Calf:	158	1.11	161
Ankle (PT):	142	1.00	144
Ankle (DP):	138	0.97	137
Digit:	107	0.75	110



Pulses

	Right		Left
Femoral:	4+		4+
Popliteal:	4+		4+
Dors. Pedis:	3+		3+
Post. Tibial:	3+		3+



Doppler






138 — **Brachial** — 142



162 — 1.14 — 169 — 1.19







158 — 1.11 — 161 — 1.13





1.00 **Ankle/Brachial Index** 1.01



Pulse Volume Recordings





138 — Brachial — 142





162 — 1.14 — 169 — 1.19



158 — 1.11 — 161 — 1.13



142 (PT) 1.00 (PT) 144 1.01
138 (DP) 0.97 (DP) 137 0.96



107 — 0.75 — TBI — 110 — 0.77

1.00 Ankle/Brachial Index 1.01

Exercise Pressures

	Rest	1	2	3	4	5	6
R Ankle (PT):	142	95	105	115	121	130	142
L Ankle (PT):	144	145	145	143	144	144	144
L Brachial:	142	158	154	150	147	145	142
R ABI:	1.00	0.60	0.68	0.77	0.82	0.90	1.00
L ABI:	1.01	0.92	0.94	0.95	0.98	0.99	1.01

Exercise Pressures

